

CAB7121

Product Information

Product SKU:	CAB7121	Gene ID:	56257	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	Refer to figures	Conjugate:	Unconjugated
Calculated MW:	74kDa	Isotype:	IgG

Immunogen Information

**Background:** Enables 7SK snRNA binding activity and RNA 5'-methyltransferase activity. Involved in RNA modification; positive regulation of protein localization to Cajal body; and positive regulation of snRNA transcription by RNA polymerase II. Located in nucleus. Part of 7SK snRNP.

Recommended Dilution:	IF/ICC, 1:50 - 1:100
Synonyms:	BCDIN3; MEPCE
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 470-689 of human MEPCE (NP_062552.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.