

LELP1 Rabbit Polyclonal Antibody



CAB7148

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

10kDa

Applications:

WB IF

Reactivity:

Human

Protein Background

Immunogen information

Gene ID:

149018

Uniprot

Q5T871

Synonyms:

LELP1

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 1-98 of human LELP1 (NP_001010857.1).

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000 IF 1:50 - 1:100

Source:

Rabbit

Isotype:

IgG


Purification:

Affinity purification

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunofluorescence analysis of A549 cells using LELP1 antibody (CAB7148). Blue: DAPI for nuclear staining.