

TBRG1 Rabbit Polyclonal Antibody



CAB7164

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

45kDa

Calculated MW:

20kDa/29kDa/44kDa

Applications:

WB

Reactivity:

Human

Protein Background

Immunogen information

Gene ID:

84897

Uniprot

Q3YBR2

Synonyms:

TBRG1; NIAM; TB-5

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 22-150 of human TBRG1 (NP_116200.2).

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG


Purification:

Affinity purification

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot analysis of extracts of various cell lines, using TBRG1 antibody (CAB7164) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (CABS014) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit (CABM00021). Exposure time: 90s.