

DZIP3 Rabbit Polyclonal Antibody



CAB7179

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

138kDa

Calculated MW:

35kDa/138kDa

Applications:

WB

Reactivity:

Human

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

9666

Uniprot

Q86Y13

Synonyms:

DZIP3; PPP1R66; UURF2; hRUL138


Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 1-230 of human DZIP3 (NP_055463.1).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot analysis of extracts of 293T cells, using DZIP3 Antibody (CAB7179) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (CABS014) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit (CABM00020). Exposure time: 15s.