

CAB7191

Product Information

| | | | | | |
|---------------------|---------|----------------------|--------|--------------------|-----------------|
| Product SKU: | CAB7191 | Gene ID: | 55081 | Size: | 20uL, 100uL |
| Clone No: | - | Host Species: | Rabbit | Reactivity: | Human,Mouse,Rat |

Additional Information

| | | | |
|-----------------------|-------|-------------------|--------------|
| Observed MW: | 49kDa | Conjugate: | Unconjugated |
| Calculated MW: | 49kDa | Isotype: | IgG |

Immunogen Information

Background: Predicted to enable DNA binding activity. Acts upstream of or within activation of cysteine-type endopeptidase activity involved in apoptotic process; apoptotic process; and regulation of apoptotic process. Predicted to be located in ciliary basal body. Predicted to be part of axoneme and intraciliary transport particle B. Predicted to be active in Golgi apparatus; centrosome; and cilium. Implicated in orofaciodigital syndrome.

Recommended Dilution: WB,1:500 - 1:2000

Synonyms: HIPPI; OFD18; MHS4R2; ESRRBL1; IFT57

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 170-429 of human IFT57 (NP_060480.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.