

CAB7206

Product Information

Product SKU:	CAB7206	Gene ID:	55825	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Mouse

Additional Information

Observed MW:	33kDa	Conjugate:	Unconjugated
Calculated MW:	33kDa	Isotype:	IgG

Immunogen Information

Background: Enables signaling receptor binding activity and trans-2-enoyl-CoA reductase (NADPH) activity. Involved in phytol metabolic process. Located in peroxisome.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	TERP; DCRRP; PVIARL; HPDHASE; SDR29C1; HSA250303; PECR
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 24-303 of human PECR (NP_060911.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.