

# ZDHHC16 Rabbit Polyclonal Antibody



CAB7226

---

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

44kDa

### Calculated MW:

34kDa/38kDa/41kDa/43kDa

### Applications:

WB

### Reactivity:

Human, Mouse, Rat

## Antibody Information

### Recommended dilutions:

WB 1:500 - 1:2000

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

## Protein Background

### Immunogen information

#### Gene ID:

84287

#### Uniprot

Q969W1

#### Synonyms:

ZDHHC16; APH2; DHHC-16

#### Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 207-296 of human ZDHHC16 (NP\_932162.1).

#### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---