
Product Information

Product SKU:	CAB7227	Gene ID:	4343	Size:	20uL, 100uL
Clone No:	ARC1817	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	114kDa	Conjugate:	Unconjugated
Calculated MW:	114kDa	Isotype:	IgG

Immunogen Information

Background: Enables 5'-3' RNA helicase activity and RNA binding activity. Involved in defense response to virus; negative regulation of transposition, RNA-mediated; and posttranscriptional regulation of gene expression. Located in P-body and cytosol. Implicated in hypertension.

Recommended Dilution: WB,1:500 - 1:1000

Synonyms: gb110; fSAP113; MOV10

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 800-900 of human MOV10 (Q9HCE1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.