

Product Information

Product SKU:	CAB7286	Gene ID:	10363	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	40kDa	Conjugate:	Unconjugated
Calculated MW:	40kDa	Isotype:	IgG

Immunogen Information

Background: Enables identical protein binding activity. Predicted to be involved in regulation of gene expression. Predicted to act upstream of or within negative regulation of neuron differentiation; negative regulation of protein sumoylation; and negative regulation of transcription by RNA polymerase II. Predicted to be active in nucleus.

Recommended Dilution:	WB,1:100 - 1:500 IF/ICC,1:50 - 1:100
Synonyms:	HMGX1; HMGXB1; HMG20A
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-347 of human HMG20A (NP_060670.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.