

Product Information

Product SKU:	CAB7296	Gene ID:	27018	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	24kDa	Conjugate:	Unconjugated
Calculated MW:	13kDa	Isotype:	IgG

Immunogen Information

Background: Enables identical protein binding activity. Predicted to be involved in signal transduction. Predicted to act upstream of or within extrinsic apoptotic signaling pathway via death domain receptors. Located in cytosol.

Recommended Dilution:	WB,1:500 - 1:2000 IF/ICC,1:50 - 1:100
Synonyms:	Bex; NADE; HGR74; NGFRAP1; DXS6984E
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-111 of human NGFRAP1 (NP_996798.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.