

CAB7313

Product Information

Product SKU:	CAB7313	Gene ID:	51499	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	12kDa	Conjugate:	Unconjugated
Calculated MW:	9kDa	Isotype:	IgG

Immunogen Information

Background: Enables p53 binding activity. Contributes to phosphatidic acid transfer activity. Involved in several processes, including DNA damage response, signal transduction by p53 class mediator; negative regulation of apoptotic process; and positive regulation of phospholipid transport. Located in mitochondrial intermembrane space and nucleoplasm. Part of protein-containing complex.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	WF-1; MDM35; P53CSV; HSPC132; TRIAP1
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-76 of human TRIAP1 (NP_057483.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.