

TRIAP1 Rabbit Polyclonal Antibody



CAB7313

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

8kDa

Applications:

WB IF

Reactivity:

Human

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000 IF 1:50 - 1:100

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

51499

Uniprot

O43715

Synonyms:

TRIAP1; HSPC132; MDM35; P53CSV; WF-1


Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 1-76 of human TRIAP1 (NP_057483.1).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunofluorescence analysis of MCF-7 cells using TRIAP1 antibody (CAB7313). Blue: DAPI for nuclear staining.