

SERPINA4 Rabbit Polyclonal Antibody



CAB7320

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

49kDa

Calculated MW:

48kDa

Applications:

WB IHC IF

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000 IHC 1:50
- 1:200 IF 1:50 - 1:100

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

5267

Uniprot

P29622

Synonyms:

SERPINA4; KAL; KLST; KST; PI-4; PI4

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 178-427 of human SERPINA4 (NP_006206.2).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images

