

**CSNK1G2 Rabbit Polyclonal Antibody**  
**CAB7326**



**Product Information**

<b>Product SKU:</b>	CAB7326	<b>Gene ID:</b>	1455	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	-	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human,Mouse,Rat

**Additional Information**

<b>Observed MW:</b>	47kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	47kDa	<b>Isotype:</b>	IgG

**Immunogen Information**

<b>Background:</b>	Enables protein serine/threonine kinase activity. Involved in peptidyl-serine phosphorylation. Located in membrane.
<b>Recommended Dilution:</b>	WB,1:500 - 1:2000
<b>Synonyms:</b>	CK1g2; CSNK1G2
<b>Purification Method:</b>	Affinity purification
<b>Immunogen:</b>	Recombinant fusion protein containing a sequence corresponding to amino acids 316-415 of human CSNK1G2 (NP_001310.3).
<b>Storage:</b>	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.