

Product Information

Product SKU:	CAB7357	Gene ID:	9912	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	110kDa	Conjugate:	Unconjugated
Calculated MW:	89kDa	Isotype:	IgG

Immunogen Information

**Background:** Enables phospholipid binding activity. Predicted to be involved in several processes, including modification of dendritic spine; negative regulation of Rac protein signal transduction; and regulation of plasma membrane bounded cell projection organization. Located in leading edge membrane.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	RICH2; NPC-A-10; ARHGAP44
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 666-772 of human ARHGAP44 (NP_055674.4).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.