

RIPK1 Rabbit Polyclonal Antibody



CAB7414

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

76kDa

Calculated MW:

70kDa/75kDa

Applications:

WB IHC IP

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000 IHC 1:50
- 1:200 IP 1:50 - 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

8737

Uniprot

Q13546

Synonyms:

RIP; RIP-1; RIP1; RIPK1

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 170-440 of human RIPK1 (NP_003795.2).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images

