

CAB7435

Product Information

Product SKU: CAB7435	Gene ID: 3104	Size: 20uL, 100uL
Clone No: -	Host Species: Rabbit	Reactivity: Human, Mouse

Additional Information

Observed MW: 76kDa	Conjugate: Unconjugated
Calculated MW: 77kDa	Isotype: IgG

Immunogen Information

Background: Enables double-stranded telomeric DNA binding activity; identical protein binding activity; and transcription cis-regulatory region binding activity. Involved in positive regulation of transcription, DNA-templated and telomere maintenance via telomere lengthening. Located in chromosome, telomeric region.

Recommended Dilution: WB, 1:100 - 1:500

Synonyms: HKR3; TZAP; ZNF855; pp9964; ZBTB48

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human ZBTB48 (NP_005332.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.