

**TM9SF1 Rabbit Polyclonal Antibody**

**CAB7461**



**Product Information**

<b>Product SKU:</b>	CAB7461	<b>Gene ID:</b>	10548	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	-	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human,Mouse

**Additional Information**

<b>Observed MW:</b>	80kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	69kDa	<b>Isotype:</b>	IgG

**Immunogen Information**

<b>Background:</b>	Predicted to be involved in protein localization to membrane. Predicted to be located in autophagosome membrane; cytoplasmic vesicle; and lysosomal membrane. Predicted to be integral component of membrane. Predicted to be active in membrane.
<b>Recommended Dilution:</b>	WB,1:500 - 1:2000
<b>Synonyms:</b>	MP70; HMP70; TM9SF1
<b>Purification Method:</b>	Affinity purification
<b>Immunogen:</b>	Recombinant fusion protein containing a sequence corresponding to amino acids 28-230 of human TM9SF1 (NP_006396.2).
<b>Storage:</b>	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.