

Product Information

Product SKU:	CAB7463	Gene ID:	23500	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	123kDa	Conjugate:	Unconjugated
Calculated MW:	123kDa	Isotype:	IgG

Immunogen Information

**Background:** Predicted to enable actin binding activity and small GTPase binding activity. Predicted to be involved in nervous system development and regulation of Wnt signaling pathway. Predicted to act upstream of or within determination of left/right symmetry. Located in extracellular exosome.

Recommended Dilution:	WB,1:500 - 1:2000 IF/ICC,1:50 - 1:200
Synonyms:	NPHS24; dJ90A20A.1; DAAM2
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-86 of human DAAM2 (NP_056160.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.