

CAB7473

Product Information

Product SKU:	CAB7473	Gene ID:	55644	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	36kDa	Conjugate:	Unconjugated
Calculated MW:	36kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable N(6)-L-threonylcarbamoyladenine synthase activity and metal ion binding activity. Involved in tRNA threonylcarbamoyladenosine modification. Located in cytoplasm; nuclear speck; and plasma membrane. Part of EKC/KEOPS complex. Implicated in Galloway-Mowat syndrome 3.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	KAE1; TCS3; GCPL1; GAMOS3; OSGEP1; PRSMG1; OSGEP
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-335 of human OSGEP (NP_060277.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.