

SAR1A Rabbit Polyclonal Antibody



CAB7476

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

26kDa

Calculated MW:

17kDa/22kDa

Applications:

WB IHC

Reactivity:

Human, Rat

Protein Background

Immunogen information

Gene ID:

56681

Uniprot

Q9NR31

Synonyms:

SAR1A; SAR1; SARA1; Sara; masra2

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 1-198 of human SAR1A (NP_064535.1).

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000 IHC 1:50
- 1:100

Source:

Rabbit

Isotype:

IgG


Purification:

Affinity purification


Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot analysis of extracts of various cell lines, using SAR1A antibody (CAB7476) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (CABS014) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit (CABM00020). Exposure time: 30s.



Immunohistochemistry of paraffin-embedded rat brain using SAR1A Antibody (CAB7476) at dilution of 1:100 (40x lens).