

Product Information

| | | | | | |
|--------------|---------|---------------|--------|-------------|-----------------|
| Product SKU: | CAB7483 | Gene ID: | 83694 | Size: | 20uL, 100uL |
| Clone No: | - | Host Species: | Rabbit | Reactivity: | Human,Mouse,Rat |

Additional Information

| | | | |
|----------------|-------|------------|--------------|
| Observed MW: | 60kDa | Conjugate: | Unconjugated |
| Calculated MW: | 60kDa | Isotype: | IgG |

Immunogen Information

Background: Predicted to enable ATP binding activity; protein serine kinase activity; and protein serine/threonine kinase activity. Predicted to be involved in protein phosphorylation. Predicted to be located in ribosome.

| | |
|-----------------------|---|
| Recommended Dilution: | WB,1:500 - 1:1000 IF/ICC,1:50 - 1:200 |
| Synonyms: | RSKL2; RPS6KL1 |
| Purification Method: | Affinity purification |
| Immunogen: | Recombinant fusion protein containing a sequence corresponding to amino acids 1-130 of human RPS6KL1 (NP_113652.2). |
| Storage: | Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3. |