

UGT2B10 Rabbit Polyclonal Antibody



CAB7570

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

65kDa/80kDa

Calculated MW:

50kDa/60kDa

Applications:

WB IHC

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000 IHC 1:50
- 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

7365

Uniprot

P36537

Synonyms:

UGT2B10; UDPGT2B10



Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 20-280 of human UGT2B10 (NP_001066.1).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot analysis of extracts of various cell lines, using UGT2B10 antibody (CAB7570) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (CABS014) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit (CABM00020). Exposure time: 90s.

Immunohistochemistry of paraffin-embedded human liver tissue using UGT2B10 antibody (CAB7570) at dilution of 1:100 (40x lens).