

UGT2B10 Rabbit Polyclonal Antibody

CAB7570



Product Information

Product SKU:	CAB7570	Gene ID:	7365	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	65kDa/80kDa	Conjugate:	Unconjugated
Calculated MW:	61kDa	Isotype:	IgG

Immunogen Information

Background:	Predicted to be involved in lipid metabolic process. Predicted to be located in endoplasmic reticulum membrane. Predicted to be integral component of membrane. Predicted to be active in intracellular membrane-bounded organelle.
Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	UDPGT2B10; UGT2B10
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 20-280 of human UGT2B10 (NP_001066.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.