

CAB7584

Product Information

Product SKU:	CAB7584	Gene ID:	22981	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Mouse

Additional Information

Observed MW:	160kDa	Conjugate:	Unconjugated
Calculated MW:	156kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable calcium ion binding activity. Predicted to be involved in microtubule anchoring at centrosome. Located in cytosol; intercellular bridge; and microtubule cytoskeleton.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	NLP; NINL
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-190 of human NINL (NP_079452.3).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.