

CAB7621

Product Information

Product SKU:	CAB7621	Gene ID:	20747	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	42kDa	Conjugate:	Unconjugated
Calculated MW:	42kDa	Isotype:	IgG

Immunogen Information

Background:	Enables RNA polymerase II-specific DNA-binding transcription factor binding activity. Acts upstream of or within several processes, including glucose homeostasis; positive regulation of apoptotic process; and regulation of transcription, DNA-templated. Located in nuclear speck. Is expressed in lung; metanephros; and pancreas. Orthologous to human SPOP (speckle type BTB/POZ protein).
Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	TEF2; Pcif1; SPOP
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-270 of mouse SPOP (NP_079563.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.