

CENPK Rabbit Polyclonal Antibody



CAB7627

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

37kDa

Calculated MW:

31kDa

Applications:

WB IF

Reactivity:

Human

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000 IF 1:50 - 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

64105

Uniprot

Q9BS16

Synonyms:

CENPK; AF5alpha; CENP-K; FKSG14; P33; Solt


Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 1-269 of human CENPK (NP_071428.2).


Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot analysis of extracts of various cell lines, using CENPK Antibody (CAB7627) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (CABS014) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit (CABM00020). Exposure time: 30s.



Immunofluorescence analysis of U2OS cells using CENPK antibody (CAB7627) at dilution of 1:100. Blue: DAPI for nuclear staining.