

Product Information

Product SKU:	CAB7656	Gene ID:	1209	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	100kDa	Conjugate:	Unconjugated
Calculated MW:	76kDa	Isotype:	IgG

Immunogen Information

**Background:** Predicted to be involved in regulation of T cell differentiation in thymus. Located in membrane.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	CLPTM1
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human CLPTM1 (NP_001285.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.