

Product Information

Product SKU:	CAB7749	Gene ID:	8395	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	71kDa	Conjugate:	Unconjugated
Calculated MW:	61kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable 1-phosphatidylinositol-4-phosphate 5-kinase activity. Predicted to be involved in regulation of phosphatidylinositol 3-kinase signaling. Predicted to act upstream of or within phosphatidylinositol biosynthetic process. Located in uropod.

Recommended Dilution:	WB,1:500 - 1:2000 IF/ICC,1:50 - 1:200
Synonyms:	MSS4; STM7; PIP5K1B
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 190-470 of human PIP5K1B (NP_003549.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.