

Product Information

|              |         |               |        |             |             |
|--------------|---------|---------------|--------|-------------|-------------|
| Product SKU: | CAB7761 | Gene ID:      | 9110   | Size:       | 20uL, 100uL |
| Clone No:    | -       | Host Species: | Rabbit | Reactivity: | Human,Mouse |

Additional Information

|                |                  |            |              |
|----------------|------------------|------------|--------------|
| Observed MW:   | Refer to figures | Conjugate: | Unconjugated |
| Calculated MW: | 133kDa           | Isotype:   | IgG          |

Immunogen Information

**Background:** Enables protein phosphatase binding activity. Involved in regulation of phosphatidylinositol dephosphorylation. Located in endosome.

|                       |   |
|-----------------------|---|
| Recommended Dilution: | IHC-P,1:50 - 1:200  |
| Synonyms:             | ZFYVE11; FYVE-DSP2; MTMR4   |
| Purification Method:  | Affinity purification   |
| Immunogen:            | Recombinant fusion protein containing a sequence corresponding to amino acids 570-820 of human MTMR4 (NP_004678.3). |
| Storage:              | Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.                  |