

CAB7806

Product Information

Product SKU:	CAB7806	Gene ID:	51647	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	18kDa	Conjugate:	Unconjugated
Calculated MW:	18kDa	Isotype:	IgG

Immunogen Information

Background: Involved in chromosome segregation; iron-sulfur cluster assembly; and protein maturation by iron-sulfur cluster transfer. Located in cytosol; nucleoplasm; and spindle. Part of CIA complex and MMXD complex.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	CIA2B; MIP18; FAM96B; CGI-128
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-163 of human FAM96B (NP_057146.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.