

ALKBH4 Rabbit Polyclonal Antibody



CAB7812

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

36kDa

Calculated MW:

4kDa/33kDa

Applications:

WB IHC

Reactivity:

Human

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000 IHC 1:50
- 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

54784

Uniprot

Q9NXW9

Synonyms:

ALKBH4; ABH4


Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 1-170 of human ALKBH4 (NP_060091.1).


Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.


Product Images




Western blot analysis of extracts of various cell lines, using ALKBH4 Antibody (CAB7812) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (CABS014) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit (CABM00020). Exposure time: 10s.



Immunohistochemistry of paraffin-embedded human stomach using ALKBH4 antibody (CAB7812) at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded human liver tissue using ALKBH4 antibody (CAB7812) at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded human lung cancer using ALKBH4 antibody (CAB7812) at dilution of 1:100 (40x lens).