

Product Information

Product SKU:	CAB7817	Gene ID:	55975	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	66kDa	Conjugate:	Unconjugated
Calculated MW:	66kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes a BTB-Kelch-related protein. The encoded protein may be involved in protein degradation. Mutations in this gene have been associated with retinitis pigmentosa 42.

Recommended Dilution:	WB,1:1000 - 1:5000
Synonyms:	CISS3; KLHL6; SBB126; PERCHING; KLHL7
Purification Method:	Affinity purification
Immunogen:	A synthetic peptide corresponding to a sequence within amino acids 200-300 of human KLHL7 (NP_061334.4).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.