

Product Information

Product SKU:	CAB7829	Gene ID:	59348	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	60kDa	Conjugate:	Unconjugated
Calculated MW:	60kDa	Isotype:	IgG

Immunogen Information

Background:	Enables DNA-binding transcription repressor activity, RNA polymerase II-specific and RNA polymerase II intronic transcription regulatory region sequence-specific DNA binding activity. Involved in negative regulation of transcription by RNA polymerase II. Located in nuclear body. Part of transcription repressor complex.
Recommended Dilution:	WB,1:500 - 1:1000
Synonyms:	ZFQR; ZBRK1; ZNF350
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 285-532 of human ZNF350 (NP_067645.3).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.