

CAB7902

Product Information

Product SKU:	CAB7902	Gene ID:	3306	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	75kDa	Conjugate:	Unconjugated
Calculated MW:	70kDa	Isotype:	IgG

Immunogen Information

Background: Enables disordered domain specific binding activity; enzyme binding activity; and unfolded protein binding activity. Involved in negative regulation of inclusion body assembly and protein refolding. Located in cytosol.

Recommended Dilution: WB,1:500 - 1:2000 IF/ICC,1:50 - 1:200

Synonyms: HSP70-2; HSP70-3; HSPA2

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 400-639 of human HSPA2 (NP_068814.2).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.