

IFNA6 Rabbit Polyclonal Antibody



CAB7903

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

22kDa

Applications:

WB IF

Reactivity:

Human

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000 IF 1:50 - 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

3443

Uniprot

P05013

Synonyms:

IFNA6; IFN-alphaK


Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 21-189 of human IFNA6 (NP_066282.1).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunofluorescence analysis of U2OS cells using IFNA6 antibody (CAB7903) at dilution of 1:100. Blue: DAPI for nuclear staining.