

Product Information

Product SKU:	CAB7948	Gene ID:	8939	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	64kDa	Conjugate:	Unconjugated
Calculated MW:	62kDa	Isotype:	IgG

Immunogen Information

Background: Enables single-stranded DNA binding activity. Involved in positive regulation of gene expression; positive regulation of transcription, DNA-templated; and transcription, DNA-templated. Acts upstream of or within positive regulation of transcription by RNA polymerase II. Located in cytoplasm and nucleus.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	FBP3; FUBP3
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 65-260 of human FUBP3 (NP_003925.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.