

Product Information

Product SKU:	CAB7978	Gene ID:	51377	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	38kDa	Conjugate:	Unconjugated
Calculated MW:	38kDa	Isotype:	IgG

Immunogen Information

Background: Enables endopeptidase inhibitor activity; proteasome binding activity; and thiol-dependent deubiquitinase. Involved in negative regulation of proteasomal ubiquitin-dependent protein catabolic process; positive regulation of smoothened signaling pathway; and protein deubiquitination. Located in cytosol; nucleolus; and nucleoplasm. Colocalizes with Ino80 complex.

Recommended Dilution: WB,1:100 - 1:500 IHC-P,1:50 - 1:200 IF/ICC,1:50 - 1:200

Synonyms: UCH37; CGI-70; INO80R; UCH-L5; UHL5

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-328 of human UHL5 (NP_001186190.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.