

Product Information

Product SKU:	CAB7992	Gene ID:	65094	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	65kDa	Conjugate:	Unconjugated
Calculated MW:	52kDa	Isotype:	IgG

Immunogen Information

Background: Enables 2-oxoglutarate-dependent dioxygenase activity. Involved in positive regulation of translational termination and protein hydroxylation. Located in cytoplasm.

Recommended Dilution:	WB,1:500 - 1:1000 IF/ICC,1:50 - 1:200
Synonyms:	JMJD4
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-463 of human JMJD4 (NP_075383.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.