

Product Information

Product SKU:	CAB7994	Gene ID:	79065	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	100kDa	Conjugate:	Unconjugated
Calculated MW:	94kDa	Isotype:	IgG

Immunogen Information

**Background:** Acts upstream of or within autophagosome assembly. Located in endosome; phagophore assembly site; and trans-Golgi network.

Recommended Dilution:	WB,1:500 - 1:1000
Synonyms:	mATG9; APG9L1; MGD3208; ATG9A
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 180-290 of human ATG9A (NP_076990.4).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.