

AHSA1 Rabbit Polyclonal Antibody



CAB8018

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

35kDa/38kDa

Calculated MW:

32kDa/38kDa

Applications:

WB

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

10598

Uniprot

O95433

Synonyms:

AHSA1; AHA1; C14orf3; hAha1; p38

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 1-338 of human AHSA1 (NP_036243.1).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
