

Product Information

Product SKU:	CAB8019	Gene ID:	23387	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	145kDa	Conjugate:	Unconjugated
Calculated MW:	145kDa	Isotype:	IgG

Immunogen Information

**Background:** Enables ATP binding activity; magnesium ion binding activity; and protein serine/threonine kinase activity. Involved in positive regulation of TORC1 signaling; positive regulation of TORC2 signaling; and protein phosphorylation. Located in nucleoplasm.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	L19; QSK; SEMDK; SIK-3; SIK3
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 60-310 of human SIK3 (NP_079440.3).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.