

Product Information

Product SKU:	CAB8022	Gene ID:	64172	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	45kDa	Conjugate:	Unconjugated
Calculated MW:	45kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable N(6)-L-threonylcarbamoyladenine synthase activity and metal ion binding activity.
Predicted to be involved in tRNA threonylcarbamoyladenosine modification. Located in mitochondrion.

Recommended Dilution:	WB,1:500 - 1:2000 IHC-P,1:50 - 1:100
Synonyms:	Qri7; OSGEPL; OSGEPL1
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-270 of human OSGEPL1 (NP_071748.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.