
Product Information

Product SKU:	CAB8030	Gene ID:	11021	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	23kDa	Conjugate:	Unconjugated
Calculated MW:	23kDa	Isotype:	IgG

Immunogen Information

Background: Enables GTPase activity; guanyl ribonucleotide binding activity; and phosphatidylinositol-4,5-bisphosphate binding activity. Involved in several processes, including endosomal transport; plasma membrane to endosome transport; and protein localization to endosome. Located in several cellular components, including clathrin-coated endocytic vesicle; clathrin-coated pit; and intercellular bridge.

Recommended Dilution: WB,1:500 - 1:1000 IHC-P,1:50 - 1:200 IF/ICC,1:50 - 1:200

Synonyms: RAY; H-ray; RAB1C; RAB35

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 100-201 of human RAB35 (NP_006852.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.