

## CAB8067

---

### Product Information

<b>Product SKU:</b> CAB8067	<b>Gene ID:</b> 55006	<b>Size:</b> 20uL, 100uL
<b>Clone No:</b> -	<b>Host Species:</b> Rabbit	<b>Reactivity:</b> Human

---

### Additional Information

<b>Observed MW:</b> 52kDa	<b>Conjugate:</b> Unconjugated
<b>Calculated MW:</b> 53kDa	<b>Isotype:</b> IgG

---

### Immunogen Information

**Background:** Enables mRNA (adenine-N1-)-methyltransferase activity; rRNA (adenine) methyltransferase activity; and tRNA (adenine-N1-)-methyltransferase activity. Involved in mRNA methylation; mitochondrial tRNA methylation; and protein homooligomerization. Located in mitochondrion.

**Recommended Dilution:** WB, 1:500 - 1:1000

**Synonyms:** TRMT61B

**Purification Method:** Affinity purification

**Immunogen:** Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human TRMT61B (NP\_060380.3).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.