

Product Information

Product SKU:	CAB8074	Gene ID:	145264	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	47kDa	Conjugate:	Unconjugated
Calculated MW:	47kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable serine-type endopeptidase inhibitor activity. Predicted to be involved in negative regulation of endopeptidase activity. Predicted to act upstream of or within negative regulation of gluconeogenesis; positive regulation of signal transduction; and regulation of lipid metabolic process. Predicted to be located in plasma membrane. Predicted to be active in extracellular space.

Recommended Dilution:	WB,1:500 - 1:2000
Synonyms:	OL-64; SERPINA12
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 20-320 of human SERPINA12 (NP_776249.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.