

RAB3IP Rabbit Polyclonal Antibody



CAB8094

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

53kDa

Calculated MW:

28kDa/33-53kDa

Applications:

WB IHC

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000 IHC 1:50
- 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

117177

Uniprot

Q96QF0

Synonyms:

RAB3IP; RABIN3; RABIN8





Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 1-150 of human RAB3IP (NP_071901.2).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot analysis of extracts of various cell lines, using RAB3IP antibody (CAB8094) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (CABS014) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit (CABM00020). Exposure time: 30s.

Immunohistochemistry of paraffin-embedded mouse brain using RAB3IP antibody (CAB8094) at dilution of 1:100 (40x lens).

Immunohistochemistry of paraffin-embedded rat brain using RAB3IP antibody (CAB8094) at dilution of 1:100 (40x lens).

Immunohistochemistry of paraffin-embedded human lung cancer using RAB3IP antibody (CAB8094) at dilution of 1:100 (40x lens).