

CAB8094

Product Information

Product SKU:	CAB8094	Gene ID:	117177	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	53kDa	Conjugate:	Unconjugated
Calculated MW:	53kDa	Isotype:	IgG

Immunogen Information

Background: Enables guanyl-nucleotide exchange factor activity and identical protein binding activity. Involved in cilium assembly; protein localization to organelle; and protein targeting to membrane. Located in centrosome; cytosol; and nucleus.

Recommended Dilution: WB,1:500 - 1:2000

Synonyms: RABIN3; RABIN8; RAB3IP

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-150 of human RAB3IP (NP_071901.2).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.