

# IFNA4 Rabbit Polyclonal Antibody



CAB8124

---

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

Refer to figures

### Calculated MW:

21kDa

### Applications:

IF

### Reactivity:

Human

## Antibody Information

### Recommended dilutions:

IF 1:50 - 1:200

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

## Protein Background

### Immunogen information

#### Gene ID:

3441

#### Uniprot

P05014

#### Synonyms:

IFNA4; IFN-alpha4a; INFA4

#### Immunogen:


Recombinant fusion protein containing a sequence corresponding to amino acids 24-189 of human IFNA4 (NP\_066546.1).

#### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---



Immunofluorescence analysis of U2OS cells using IFNA4 antibody (CAB8124) at dilution of 1:100. Blue: DAPI for nuclear staining.