

Product Information

Product SKU:	CAB8143	Gene ID:	5540	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	42kDa	Conjugate:	Unconjugated
Calculated MW:	42kDa	Isotype:	IgG

Immunogen Information

Background: Enables pancreatic polypeptide receptor activity and peptide hormone binding activity. Involved in G protein-coupled receptor signaling pathway. Located in membrane.

Recommended Dilution:	WB, 1:500 - 1:2000
Synonyms:	Y4; PP1; PPYR1; NPY4-R; NPY4R
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-80 of human NPY4R (NP_005963.4).
Storage:	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.